English – INSTRUCTION GUIDE

EYEPRIM™ - Conjunctival impression device

MADE IN FRANCE

OPIA Technologies SAS
17, rue Moreau 75012
Paris, France
Email : sales@opiatech.com

REF OS1


55°C


STERILE EO

Sterility guaranteed unless unit pack has been damaged or opened. Do not reuse. Do not re-sterilize. Discard after use.

COMPOSITION

EYEPRIM™ Device

Membrane position on one of the 4 quadrants of the eyeball

EYEPRIM™ consists of 3 main elements:

- A rigid body, with handles, allowing a good handling.
- A flexible membrane inserted at the distal end of the device and protected within it, to be applied to the conjunctiva.
- A push button at the other end, to bring the membrane in contact with the eye and eject the membrane.

INDICATIONS

For patients with ocular surface disease whose diagnosis (dry eye, allergy, infection...) can be clarified by molecular biology analysis. EYEPRIM™ allows the collection of a quadrant on the patient's eye: the choice of quadrant is left to the physician choice.

CONTRAINDICATIONS


For children aged less than 3 years old, it is strongly recommended to check that the applicator fits with the ocular globe's size.

INSTRUCTIONS FOR USE

1 - The patient can be indifferently in a vertical (head slightly leaning behind) or horizontal position during the sampling act. The patient must look up or down according to the selected quadrant. For sensitive patients, it may be wise to practice a topical anaesthesia of the eye to be treated. Notice that anaesthesia could reduce the sampling performances. In case of PCR analysis, do not use anaesthetic or other topical substances that could inhibit PCR efficiency.

2 - Open the package.

3 - Keep the eyelids open with the thumb and forefinger, or with the help of a blepharostat.

4 – Position the distal end of the device upon the eye conjunctiva of the patient. The inner radius can be positioned relative to the limbus.

5 – Press on the push-button until reaching the stop to bring the membrane into contact with the conjunctiva. Hold from 2 to 3 seconds.

6 – **Release** the pressure before removing the device from the eye to protect the sample from the environment.

* Notice that the membrane is stretched onto its support

7 – Place the distal end of the device over the storage container. The inner diameter of the storage container should be at least 2 cm.

8 – Eject the membrane by pushing hard onto the push-button until depression of the latter.

9 - Discard the device according to hospital standard operating procedures.

- . In case where the membrane remains on the patient's eye, remove it gently with forceps by grabbing the membrane's tab.
- . Medical staff only.
- . Handle with gloves.
- . To monitor the evolution of pathology, it is recommended to always take the same quadrant.
- . The membrane cannot be made transparent.
- . The device is not intended to be reused: this could cause cross infection.

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